



Komet Medical

Pins für die Knieendoprothetik.

Pins for knee joint endoprosthesis.

Mit einer neuen Produktlinie für die Knieendoprothetik setzt Komet Medical einmal mehr unter Beweis, Vollsortimenter im Bereich motorbetriebener Werkzeuge für chirurgische Antriebssysteme zu sein. Pins spielen eine elementare Rolle bei der Fixierung der Knieschablone.

Das Komet Medical Sortiment umfasst verschiedene Arten von Pins und Screws, die in steriler und unsteriler Verpackung ab Lager verfügbar sind.

Once again, our new product range for knee joint endoprosthesis has shown that Komet Medical is a full-range supplier of motor-driven tools for surgical power systems. Pins play a decisive role in the fixation of knee templates. The Komet Medical range comprises different types of pins, all of which are available from stock in sterile or non-sterile packaging.



Guide Pins | Guide Screws



Bohrpin / Pin	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GP003.070.023 STE	2,3	70,0	-	Rundschaft / Roundshank ●	-	-
GP003.070.023	2,3	70,0	-	Rundschaft / Roundshank ●	-	-




Bohrpin / Pin	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GP002.089.032 STE	3,2	89,5	-	4-kant / Square ■	2,25	GS900.100.010,
GP002.089.032	3,2	89,5	-	4-kant / Square ■	2,25	GS900A.040.000



Universal Pin	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GP001.089.032 STE	3,2	89,0	70,0	4-kant / Square ■	2,25	GS900.100.010,
GP001.089.032	3,2	89,0	70,0	4-kant / Square ■	2,25	GS900A.040.000



Universal Pin 	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GP001A.089.032 STE	3,2	89,0	70,0	4-kant / Square ■	2,25	GS900.100.010,
GP001A.089.032	3,2	89,0	70,0	4-kant / Square ■	2,25	GS900A.040.000



Schraubpin / Screw 	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS905.030.032 STE	3,2	30,0	26,0	6-kant / Hexagon ●	3,5	-
GS905.030.032	3,2	30,0	26,0	6-kant / Hexagon ●	3,5	-



Schraubpin / Screw 	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS901.033.032 STE	3,2	33,0	30,0	6-kant / Hexagon ●	2,5	GS900.100.000,
GS901.033.032	3,2	33,0	30,0	6-kant / Hexagon ●	2,5	GS900.130.000



Schraubpin / Screw 	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS901.050.032 STE	3,2	50,0	46,0	6-kant / Hexagon ●	2,5	GS900.100.000,
GS901.050.032	3,2	50,0	46,0	6-kant / Hexagon ●	2,5	GS900.130.000



Schraubpin / Screw 	Ø (mm)	L (mm)	GL / TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS901A.050.032 STE	3,2	50,0	46,0	6-kant / Hexagon ●	2,5	GS900.100.000,
GS901A.050.032	3,2	50,0	46,0	6-kant / Hexagon ●	2,5	GS900.130.000



Schlüsselweite
Wrench width



4-kant
Square



6-kant
Hexagon



Rundschaft
Roundshank



magnetisierbarer Stahl
magnetizable steel



Schraubpin / Screw		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS905.050.032 STE		3,2	50,0	46,0	6-kant/Hexagon	3,5	-
GS905.050.032		3,2	50,0	46,0	6-kant/Hexagon	3,5	-



Schraubpin / Screw		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS903.062.032 STE		3,2	62,5	46,0	4-kant/Square	2,25	GS900.100.010,
GS903.062.032		3,2	62,5	46,0	4-kant/Square	2,25	GS900A.040.000



Schraubpin / Screw		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS903A.062.032 STE		3,2	62,5	46,0	4-kant/Square	2,25	GS900.100.010,
GS903A.062.032		3,2	62,5	46,0	4-kant/Square	2,25	GS900A.040.000



Pin (lang) / Long Pin		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS902.150.032 STE		3,2	150,0	51,0	Rundschaft/Roundshank	-	-
GS902.150.032		3,2	150,0	51,0	Rundschaft/Roundshank	-	-



Pin (lang) / Long Pin		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS904.100.030 STE		3,0	100,0	26,0	4-kant/Square	2,25	GS900.100.010,
GS904.100.030		3,0	100,0	26,0	4-kant/Square	2,25	GS900A.040.000




Pin (lang) / Long Pin		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS904.150.030 STE		3,0	150,0	26,0	4-kant/Square	2,25	GS900.100.010,
GS904.150.030		3,0	150,0	26,0	4-kant/Square	2,25	GS900A.040.000




Pin (lang) / Long Pin		Ø (mm)	L (mm)	GL/TL (mm)	Schaft / Shank	Schlüsselweite / Wrench width (mm)	Kompatibel mit Adapter / Compatible with adapter
GS904.100.040 STE		4,0	100,0	25,0	4-kant/Square	2,25	GS900.100.010,
GS904.100.040		4,0	100,0	25,0	4-kant/Square	2,25	GS900A.040.000




Pin Adapter	L (mm)	Schlüsselform/ Wrench shape	Schlüsselweite / Wrench width (mm)	Schaft / Shank
GS900.100.000	100,0	6-kant/Hexagon 	2,5	Zimmer/Hudson




Pin Adapter	L (mm)	Schlüsselform/ Wrench shape	Schlüsselweite / Wrench width (mm)	Schaft / Shank
GS900.100.010	100,0	4-kant/Square 	2,25	Zimmer/Hudson



Pin Adapter	L (mm)	Schlüsselform/ Wrench shape	Schlüsselweite / Wrench width (mm)	Schaft / Shank
GS900.130.000	130,0	6-kant/Hexagon 	2,5	Zimmer/Hudson



Pin Adapter	L (mm)	Schlüsselform/ Wrench shape	Schlüsselweite / Wrench width (mm)	Schaft / Shank
GS900A.040.000	40,0	4-kant/Square 	2,25	AO